


## 4.5mm Distal Femur Locking Plate





# PERI-LOC Periarticular Locked Plating System

## Distal Femur Plate Surgical Technique

### Contents

Product Overview .....	2
Indications .....	2
Design Features .....	3
Patient Positioning .....	4
Incision .....	4
4.5mm Distal Femur Locking Plate Surgical Technique	
Plate Selection .....	6
Articular Reduction and Provisional Fixation .....	7
Plate Positioning .....	7
Screw Insertion .....	9
Catalog Information .....	13

### Nota Bene

The technique description herein is made available to the healthcare professional to illustrate the author's suggested treatment for the uncomplicated procedure. In the final analysis, the preferred treatment is that which addresses the needs of the specific patient.

# Product Overview

The PERI-LOC Periarticular Locked Plating System from Smith & Nephew, Inc. offers the advantages of locked plating with the flexibility and benefits of traditional plating in one system. Utilizing both locking and non-locking screws, PERI-LOC offers a construct that resists angular (e.g. varus/valgus) collapse while simultaneously acting as an effective aid to fracture reduction. A simple and straightforward instrument set features one screwdriver, standardized drill bits, and color-coded instrumentation, making PERI-LOC efficient and easy to use.

All PERI-LOC implants are manufactured using the highest quality 316L stainless steel for strength and durability.

The anatomical bow and precontour of the 4.5mm Distal Femur Locking Plate provides an excellent fit against the surface of the bone.

Condylar scallops on the distal end of the plate allow easy placement of lag screws outside the plate for fixation of articular fractures.

Each screw hole will accept one of four different screws allowing you to customize the screw configuration depending on the individual needs of the fracture.

- 4.5mm Self-Tapping Cortex Screw (Non-Locking)
- 4.5mm Locking Self-Tapping Cortex Screws
- 5.7mm Cannulated Locking Screw
- 6.5mm Partially Threaded Cancellous Screw




## Indications

The PERI-LOC Periarticular Locked Plating System is used for adult and pediatric patients as indicated for pelvic, small, and long bone fracture fixation. Indications for use include fractures of the tibia, fibula, femoral condyle, pelvis, acetabulum, metacarpals, metatarsals, humerus, ulna, middle hand and middle foot bones (particularly in osteopenic bone); treatment of the calcaneus; hip arthrodesis, and provisional hole fixation.

Components in the PERI-LOC Periarticular Locked Plating System are for single use only.

# Design Features



**Each of the holes can accept one of four different screws:**



4.5mm Self-Tapping Cortex Screw (Non-Locking)  
Cat. No. 7182-60XX



4.5mm Locking Self-Tapping Cortex Screw  
Cat. No. 7182-70XX



5.7mm Cannulated Locking Screw  
Cat. No. 7182-80XX



6.5mm Partially Threaded Cancellous Screw  
Cat. No. 7182-81XX

All screws use 3.5mm Hexdriver.


# Patient Positioning

Place the patient in a supine position on a radiolucent table. A small bump can be used under the ipsilateral hip. The entire leg and lateral hip region should be prepped and draped to allow proximal extension of the surgical exposure if necessary. A sterile tourniquet can be used, especially for distal fractures. Confirm that an unhindered lateral and AP view under fluoroscopy can be obtained.


Obtain gross metaphyseal alignment using manual traction or skeletal distraction.

## Incision


The incision illustrated below is indicated for the following fractures:



Articular simple, metaphyseal simple (33-C1) (optional)




Articular simple, metaphyseal multifragmentary (33-C2)




Multifragmentary articular fracture (33-C3)


OTA Fracture Classification courtesy of the Orthopaedic Trauma Association. For more information go to [www.ota.org](http://www.ota.org)




The incision illustrated below is indicated for the following fractures:



Extra articular (33-A)




Partial articular (33-B)




Articular simple, metaphyseal simple (33-C1) (optional)

OTA Fracture Classification courtesy of the Orthopaedic Trauma Association. For more information go to [www.ota.org](http://www.ota.org)





## Exposure



# 4.5mm Distal Femur Locking Plate Surgical Technique

## Plate Selection

Using the PERI-LOC 4.5mm Distal Femur Locking Plate Preoperative Template, determine the appropriate length plate for the fracture. In general, a longer plate allows for better mechanical advantage over a shorter plate. It is recommended that when selecting plate length, allow for five screw holes above the most proximal aspect of the fracture.




PERI-LOC 4.5mm Distal Femur Plating Preoperative Template  
Cat. No. 7118-0915



## Articular Reduction and Provisional Fixation


It is important that articular fracture reduction be obtained prior to placement of locking screws. Temporarily secure articular fragments by using K-Wires and/or Reduction Forceps. Place provisional and/or definitive fixation peripheral to the condylar contour of the plate. **3.5mm and 4.5mm Self-Tapping Cortex Screws (Non-Locking) can be nested peripherally in the contours of the plate.**



**NOTE:** If a **posterior Hoffa fracture** is present, fixation can be obtained by placing 3.5mm cortex screws or 4.0mm cancellous screws from anterior to posterior. Be sure to countersink the screw heads by using the Large Fragment Countersink so that the screw heads rest below the level of articular cartilage.

## Position Plate

Position the PERI-LOC 4.5mm Distal Femur Locking Plate by matching the contour of the plate to the distal portion of the lateral femur. Place the long (metaphyseal) Provisional Fixation Pin through the center hole of the distal cluster.



Reduction Forceps  
Cat. No. 7117-0044



Large Fragment  
Countersink  
Cat. No. 7117-3353



3.5mm Provisional  
Fixation Pin 40mm  
Cat. No. 7117-3325




3.5mm Self-Tapping  
Cortex Screw  
(Non-Locking)  
Cat. No. 7182-40XX



4.5mm Self-Tapping  
Cortex Screw  
(Non-Locking)  
Cat. No. 7182-60XX

Attach the 4.5mm/5.7mm Locking Screw Guide to any of the distal holes and insert the 2.0mm K-Wire Locking Guide Insert (blue) which accepts the 2.0mm K-Wire (guide wire). This K-Wire can be re-directed if necessary until it is parallel to the joint in the AP view. Loosening of the Provisional Fixation Pin may be necessary to redirect the K-Wire parallel to the joint.



For correct coronal alignment the K-Wire (guide wire) must be parallel to the joint in the AP view.

Advance the K-Wire until it reaches the medial wall of the femoral condyle. Measure for screw length by placing the 5.7mm Cannulated Depth Gauge against the end of the 2.0mm K-Wire Locking Guide Insert for proper measurement.



4.5mm/5.7mm  
Locking Screw Guide  
Cat. No. 7117-3539



2.0mm K-Wire  
Locking Guide Insert  
Cat. No. 7117-3531



2.0mm X 228mm K-Wire  
Cat. No. 7117-3361




5.7mm Cannulated  
Depth Gauge  
Cat. No. 7117-3526

## Screw Insertion


Remove the Locking Screw Guide with Insert and insert the appropriate length 5.7mm Cannulated Locking Screw over the K-Wire and into the bone using the 3.5mm Cannulated Hexdriver Shaft.

**Note:** The 5.7mm Cannulated Screws are self-drilling and self-tapping, making predrilling unnecessary in most cases. However, if predrilling is necessary, drill the near cortex using the 4.5mm Cannulated Drill Bit with Quick Connect.



Obtain sagittal alignment and confirm with a lateral radiograph. Reference the position of the plate to Blumensaat's Line and the subchondral margin of the trochlear groove.

Center the plate on the lateral aspect of the femur and apply a short (diaphyseal) Provisional Fixation Pin in the second most proximal hole.



5.7mm Cannulated  
Locking Screw  
Cat. No. 7182-80XX



3.5mm Cannulated  
Hexdriver Shaft  
Cat. No. 7117-3536




4.5mm Cannulated Drill  
Bit with Quick Connect  
Cat. No. 7117-3508




3.5mm x 40mm  
Provisional Fixation Pin  
Cat. No. 7117-3325

Center the plate on the distal diaphyseal fracture fragment and provisionally fix the plate close to the fracture. Obtain final confirmation of fracture alignment and implant position.



The remaining condylar screws can be either 5.7mm Cannulated Locking Screws or 4.5mm Locking Self-Tapping Cortex Screws. To implant 4.5mm Locking Self-Tapping Cortex Screws, predrill with the 3.5mm Drill Bit with Quick Connect through the 4.5mm/5.7mm Locking Screw Guide and 3.5mm Locking Drill Guide Insert (red), stopping short of the medial cortex.

**Note:** Locking screws can be inserted using a powered drill system but should be tightened by hand. Tightening screws using a powered drill system may cause loss of reduction or expose the screw heads to excess torque.



4.5mm Locking  
Screw  
Cat. No. 7182-70XX



3.5mm Drill Bit with  
Quick Connect  
Cat. No. 7117-3505



4.5mm/5.7mm Locking  
Screw Guide  
Cat. No. 7117-3539



3.5mm Locking Drill  
Guide Insert (Red Insert)  
Cat. No. 7117-3530

The use of at least one 5.7mm Cannulated Locking Screw is recommended in the distal fragment.

**Note:** It may be necessary to use a unicondylar screw in the most distal hole to avoid joint impingement.



Proceed with definitive fixation of the shaft and the condylar portions with appropriate screw selections. If a combination of non-locking screws and locking screws is necessary, insert the non-locking cortex screws first to reduce the plate to bone or to compress the fracture before locking screws are inserted.

Pre-drill for the self-tapping cortex screws (non-locking) using the 3.5mm Drill Bit through the 3.5mm Compression or Neutral Locking Hole Insert (green or gold round drill guide inserts). Measure for length using the calibrations on the 3.5mm Drill Bit or Large Fragment Screw Depth Gauge and insert the appropriate length 4.5mm Self-Tapping Cortex Screw (Non-Locking) using the 3.5mm Hexdriver with Quick Connect.



3.5mm Neutral Locking Hole Insert  
Cat. No. 7117-3521



3.5mm Compression Locking Hole Insert  
Cat. No. 7117-3522



Large Fragment Screw Depth Gauge  
Cat. No. 7117-3331



4.5mm Self-Tapping Cortex Screw  
Cat. No. 7182-60XX



3.5mm Hexdriver Shaft with Quick Connect  
Cat. No. 7117-3537


After compression is obtained, thread the 4.5mm/5.7mm Locking Screw Guide with 3.5mm Locking Drill Guide Insert (red) in screw holes where 4.5mm Locking Self-Tapping Cortex Screws are needed. Again, drill using the 3.5mm Drill Bit and read measurement from drill bit or use the Large Fragment Screw Depth Gauge. Insert the appropriate length 4.5mm Locking Self-Tapping Cortex Screw using the 3.5mm Hexdriver with Quick Connect. Make sure all screws are tight before closing the wound.



Final lateral view



Final AP view



# Catalog Information – 4.5mm Distal Femur Plates

## 4.5mm Distal Femur Locking Plates

Cat. No.	Length	Quantity in Set
7182-0006	6H Left 155mm	1
7182-0008	8H Left 193mm	1
7182-0010	10H Left 230mm	1
7182-0013	13H Left 286mm	1
7182-0016	16H Left 342mm	1
7180-0019	19H Left 399mm	0
7182-0106	6H Right 155mm	1
7182-0108	8H Right 193mm	1
7182-0110	10H Right 230mm	1
7182-0113	13H Right 286mm	1
7182-0116	16H Right 342mm	1
7180-0119	19H Right 399mm	0



## Small Outer Case – 2.4”

Cat. No. 7112-9401

## Lid for Outer Cases

Cat. No. 7112-9402

## Plate Tray

Cat. No. 7117-0323

# Catalog Information – Large Fragment System Screws

## Large Fragment System 4.5mm Self-Tapping Cortex Screws (Non-Locking)



Cat. No.	Length	Quantity in Set
7182-6014	14mm	4
7182-6016	16mm	4
7182-6018	18mm	4
7182-6020	20mm	6
7182-6022	22mm	6
7182-6024	24mm	6
7182-6026	26mm	6
7182-6028	28mm	6
7182-6030	30mm	10
7182-6032	32mm	10
7182-6034	34mm	10
7182-6036	36mm	10
7182-6038	38mm	10
7182-6040	40mm	10
7182-6042	42mm	6
7182-6044	44mm	4
7182-6046	46mm	4
7182-6048	48mm	4
7182-6050	50mm	4
7182-6052	52mm	4
7182-6054	54mm	4
7182-6056	56mm	4
7182-6058	58mm	4
7182-6060	60mm	4
7182-6062	62mm	4
7182-6064	64mm	4
7182-6066	66mm	4
7182-6068	68mm	4
7182-6070	70mm	4
7182-6072	72mm	4
7182-6074	74mm	4
7182-6076	76mm	4
7182-6078	78mm	4
7182-6080	80mm	4
7182-6085	85mm	4
7182-6090	90mm	2
7182-6095	95mm	2
7182-6100	100mm	2
7180-6105	105mm	0
7180-6110	110mm	0
7180-6115	115mm	0
7180-6120	120mm	0
7180-6125	125mm	0
7180-6130	130mm	0



## Large Fragment System 4.5mm Locking Self-Tapping Cortex Screws



Cat. No.	Length	Quantity in Set
7182-7010	10mm (Blunt Tip)	4
7182-7012	12mm (Blunt Tip)	4
7182-7014	14mm	4
7182-7016	16mm	4
7182-7018	18mm	4
7182-7020	20mm	6
7182-7022	22mm	6
7182-7024	24mm	6
7182-7026	26mm	6
7182-7028	28mm	6
7182-7030	30mm	10
7182-7032	32mm	10
7182-7034	34mm	10
7182-7036	36mm	10
7182-7038	38mm	10
7182-7040	40mm	10
7182-7042	42mm	6
7182-7044	44mm	4
7182-7046	46mm	4
7182-7048	48mm	4
7182-7050	50mm	4
7182-7052	52mm	4
7182-7054	54mm	4
7182-7056	56mm	4
7182-7058	58mm	4
7182-7060	60mm	4
7182-7062	62mm	4
7182-7064	64mm	4
7182-7066	66mm	4
7182-7068	68mm	4
7182-7070	70mm	4
7182-7072	72mm	4
7182-7074	74mm	4
7182-7076	76mm	4
7182-7078	78mm	4
7182-7080	80mm	4
7182-7085	85mm	4
7182-7090	90mm	2
7182-7095	95mm	2
7182-7100	100mm	2
7180-7105	105mm	0
7180-7110	110mm	0
7180-7115	115mm	0
7180-7120	120mm	0
7180-7125	125mm	0
7180-7130	130mm	0

## Large Fragment System 5.7mm Cannulated Locking Screws



Cat. No.	Length	Quantity in Set
7182-8020	20mm	3
7182-8025	25mm	3
7182-8030	30mm	3
7182-8035	35mm	3
7182-8040	40mm	3
7182-8045	45mm	3
7182-8050	50mm	3
7182-8055	55mm	5
7182-8060	60mm	5
7182-8065	65mm	5
7182-8070	70mm	5
7182-8075	75mm	5
7182-8080	80mm	5
7182-8085	85mm	3
7182-8090	90mm	3
7182-8095	95mm	3
7182-8100	100mm	3
7180-8105	105mm	0
7180-8110	110mm	0
7180-8115	115mm	0
7180-8120	120mm	0

## 6.5mm Partially Threaded Cancellous Screws



Cat. No.	Length	Quantity in Set
7182-8150	50mm	4
7182-8155	55mm	4
7182-8160	60mm	4
7182-8165	65mm	4
7182-8170	70mm	4
7182-8175	75mm	4
7182-8180	80mm	4
7182-8185	85mm	4
7182-8190	90mm	4
7182-8195	95mm	4
7182-8200	100mm	4
7180-8205	105mm	0
7180-8210	110mm	0

## Washers

Cat. No.	Length	Quantity in Set
7114-3110	10mm O.D.	6
7114-3113	13mm O.D.	6



# Catalog Information – Large Fragment System Instruments

## Sharp Hook

Cat. No. 7117-0043



## Wire Bending Pliers, 140mm Length

Cat. No. 7117-0063



## Large Fragment Screw Depth Gauge

Cat.No. 7117-3331



## 5.7mm Cannulated Depth Gauge

Cat.No. 7117-3526



## Large Fragment Countersink

Cat.No. 7117-3353



## Universal Plate Bending Irons

Cat.No. 7117-3367



## Hohmann Retractor Long, 15mm Width

Cat.No. 7117-3393



## Universal Drill Guide Handle

Cat.No. 7117-3349



## 2.0mm Wire/Drill Insert

Cat.No. 7117-3517



## 3.5mm Drill Guide Insert

Cat.No. 7117-3513



## 2.0mm Parallel Wire/Drill Guide

Cat.No. 7117-3516



## 4.5mm Drill Guide Insert

Cat.No. 7117-3520



## 3.5mm Neutral Locking Hole Insert

Cat.No. 7117-3521



## 3.5mm Compression Locking Hole Insert

Cat.No. 7117-3522



## 3.5mm Neutral Slot Insert

Cat.No. 7117-3519



## 3.5mm Compression Slot Insert

Cat.No. 7117-3518



## 4.7mm Hexdriver

Cat.No. 7117-3540



Cannulated Bending Irons for K-Wires  
Cat.No. 7117-3527



Cannulated AO to Trinkle Adaptor  
Cat.No. 7117-3528



4.5/5.7mm Locking Screw Guide  
Cat.No. 7117-3539



2.0mm K-Wire Locking Guide Insert  
Cat.No. 7117-3531



3.5mm Locking Drill Guide Insert  
Cat.No. 7117-3530



4.5mm Locking Drill Guide Insert  
Cat.No. 7117-3532



3.5mm Locking Drill Guide – One Piece  
Optional  
Cat. No. 7117-3451



4.5mm Locking Drill Guide – One Piece  
Optional  
Cat. No. 7117-3541



Large Fragment Guide Removal Assembly  
Cat.No. 7117-3550



Large Screwdriver Handle  
Cat.No. 7117-3547



Tear Drop Handle Screwdriver with Quick Connect  
Cat.No. 7117-3543



Small T-Handle, Quick Coupling  
Cat.No. 7117-3542



3.5mm Hexdriver Shaft  
with AO Quick Connect  
Cat.No. 7117-3537



3.5mm Cannulated Hexdriver Shaft  
Cat.No. 7117-3536



# Catalog Information – Large Fragment System Forceps Tray Instruments

## Self Centering Reverse Verbrugge

Cat. No.	Description
7117-3544	190mm
7117-3545	240mm
7117-3546	280mm



## Reduction Forceps with Ratchet, 205mm

Cat. No. 7117-0044



## Reduction Forceps with Speed Knob, 240mm

Cat. No. 7117-0050



## Socket Wrench with Universal Joint

Cat. No. 7117-0143



## Articulated Tension Device with Gauge

Cat. No. 7117-0145



## Lamina Spreader

Cat. No. 7117-3365



## Reduction Forceps with Ratchet-Bowed, 205mm

Cat. No. 7117-3370



## Reduction Forceps with Ratchet, 240mm

Cat. No. 7117-3371



## Reduction Forceps with Points, Broad

Cat. No. 7117-3377



## Reduction Forceps with Serrated Jaw

Cat. No. 7117-3378



# Catalog Information – Large Fragment System Trays

## PERI-LOC Large Fragment Instrument Tray

Cat.No. 7117-0327

## Small Outer Case – 2.4”

Cat. No. 7112-9401

## Lid for Outer Cases

Cat. No. 7112-9402

## PERI-LOC Forceps Tray

Cat. No. 7117-0326

# Catalog Information – Large Fragment System Disposables

## K-Wires with Trocar Point and Threaded Pins

Cat. No.	Description	Quantity in Set
7116-1020	2.0mm x 150mm	6
7117-3361	2.0mm x 228mm	6



## Taps with Quick Connect

Cat. No.	Description	Quantity in Set
7117-3319	4.5mm	2
7117-3509	6.5mm Cancellous	2



## Provisional Fixation Pins

Cat. No.	Description	Quantity in Set
7117-3324	3.5mm x 18mm	4
7117-3325	3.5mm x 40mm	4



## Drill Bits with Quick Connect

Cat. No.	Description	Quantity in Set
7117-3504	3.5mm Short	2
7117-3505	3.5mm	2
7117-3506	4.5mm	2
7117-3507	4.5mm Short	2
7117-3508	4.5mm Cannulated	2





**Orthopaedics**

Smith & Nephew, Inc.  
1450 Brooks Road  
Memphis, TN 38116  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)

Telephone: 901-396-2121  
Information: 1-800-821-5700  
Orders/inquiries: 1-800-238-7538

The following statement is required by the U.S. FDA: **WARNING:** This device is not approved for screw attachment or screw fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.