

Designed with patient wellbeing in mind

ALLEVYN[®] LIFE dressing was designed with the input of both patients and clinicians from around the world to ensure a product that delivers the best wound care experience.

- The quadrilobe shape combined with the wider border helps ensure ALLEVYN LIFE fits securely to the body even during a shower⁷
- The protective masking layer minimises the visual impact of exudate strikethrough (*in vitro*) for patient confidence⁹
- The lock-away layer minimises the risk of fluid leakage, providing patient peace of mind⁸
- The cushioning layer helps maximise protection from every day knocks and bumps¹⁰
- The gentle silicone adhesive makes the dressing easy to apply, remove and reposition as well as helping ensure patient comfort during wear and on removal^{7,11}

Whilst maintaining patient satisfaction:

Patient satisfaction results from a recent clinical evaluation (N=148)⁷

- 71% were satisfied with odour control
- 70% dressing retention
- 69% leakage protection
- 84% comfort during wear
- 66% exudate masking
- 81% pain on dressing removal

ALLEVYN LIFE is also available in specially designed Sacrum and Heel shapes



With ALLEVYN[®] LIFE, you can treat wounds efficiently and effectively

Designed to save nurses time by...

- ✓ Reducing the need for frequent dressing changes⁴⁻⁶
- ✓ Reducing visit frequency⁴

...whilst maintaining patient satisfaction⁴

Try ALLEVYN LIFE in your practice.



ALLEVYN LIFE		
S&N Code	Dressing Size (cm)	Carton
66801067	10.3cm x 10.3cm	10
66801068	12.9cm x 12.9cm	10
66801069	15.4cm x 15.4cm	10
66801070	21cm x 21cm	10
66801304	Heel 25cm x 25.2cm	5
66801306	Small Sacrum 17.2cm x 17.5cm	10
66801307	Large Sacrum 21.6cm x 23cm	10



Visit www.allevynlife.com
or scan the QR code

Supporting healthcare professionals for over 150 years

Wound Management
Smith & Nephew
Medical Ltd
101 Hessle Road
Hull HU3 2BN
UK


T +44 (0) 1482 225181
F +44 (0) 1482 328326

[®]Trademark of Smith & Nephew
All Trademarks acknowledged
©May 2016 Smith & Nephew
69354

References

1. Royal College of Nursing. Overstretched. Under-resourced. The UK nursing labour market review 2012. RCN: London, UK, 2012.
2. O'Keeffe M. Evaluation of a community-based wound care programme in an urban area. Poster presented at EWMA Conference 2006.
3. MedMarket Diligence, LLC. Wound prevalence and wound management, 2012-2020. Available at: <http://blog.medilgence.com/2013/01/29/wound-prevalence-and-wound-management-2012-2020/> [Accessed February 2016].
4. Joy H *et al.* A collaborative project to enhance efficiency through dressing change practice. *Journal of Wound Care* 2015;24(7):312,314-7.
5. Stephen-Haynes J *et al.* An appraisal of the clinical performance and economic benefits of a silicone foam in a large UK primary care organisation. *Journal of Community Nursing* 2013; 6.
6. Simon D and Bielby A. A structured collaborative approach to appraise the clinical performance of a new product. *Wounds UK* 2014;10(3):80-87.
7. Rossington A *et al.* Clinical performance and positive impact on patient wellbeing of ALLEVYN LIFE. *Wounds UK* 2013;9(4):91-95.
8. Smith & Nephew Data on file report: DS/12/130.DOF. Simulated wound model testing of ALLEVYN LIFE and Mepilex[™] Border.
9. Smith & Nephew Data on file report: DS/12/129/ DOF. Simulated wound model testing of ALLEVYN LIFE and Mepilex[™] Border.
10. Smith & Nephew Data on file report: DS/12/185/DOF. Impact protection properties of ALLEVYN LIFE, Mepilex[™] Border and Biatain[™] Silicone.
11. Smith & Nephew Data on file report: OR DOF020. Performance of silicone dressings – ALLEVYN LIFE and Mepilex[™] Border.

Care takes time. Now you have it.



smith&nephew
ALLEVYN[®] LIFE

Designed for people who
happen to be patients

Nurses today are under more pressure than ever...



...due to factors such as:¹



and research shows that...



With the number of wounds requiring treatment expected to rise in the coming years,³ **it's time to take a different approach to dressing changes.**

"I barely have enough time to care for all my patients as it is, and that's before I worry about wound care."

Nurse

* Based on staffing levels at the time of the study

ALLEVYN[®] LIFE – Saving nurses time through fewer dressing changes...

ALLEVYN LIFE is an innovative, multi-layered foam dressing that has been shown to help save time by reducing the number of dressing changes.⁴⁻⁶

The unique five layer design provides optimal fluid handling:⁷

- Shower-proof, highly breathable film
- Built-in discretion (masking layer)
- Highly absorbent, lock-away core
- Hydrocellular foam
- Silicone adhesive wound contact layer



The hyper-absorbent lock-away layer helps ensure no leakage of fluid.⁸



Percentage of dressings with leakage on a 7 day wound model test (*in-vitro*) of ALLEVYN LIFE vs. Mepilex[™] Border (N=22)⁸

The unique change indicator has been shown to increase the time between dressing changes.⁵



...and reduced visit frequency

In a recent evaluation, ALLEVYN[®] LIFE significantly reduced the number of dressing changes through increased mean wear time.^{6†}

Foam Dressings



* when used with appropriate training/education and effective promotion of practice change.

Reducing dressing change visits

Clinical studies show that with ALLEVYN LIFE you can treat wounds efficiently and effectively and reduce dressing change visit frequency.^{4,5}

Evaluations of ALLEVYN LIFE have shown:

- 50% reduction in the number of visits for dressing changes⁴
 - 2,500 potential nurse working days could be saved when treating patients needing frequent dressing changes (≥ 3 visits/week; N=469) over a period of one year^{5†}
- [†] estimated from wear time, trust demographic and audit data.

98% of clinicians said they would recommend ALLEVYN LIFE for use in their practice.⁶