

Arthroscopes and VideoArthroscopes



Excellent engineering, efficient design, and exceptional performance—the qualities you expect from the leader in Endoscopic imaging.



The most technically advanced product portfolio for Arthroscopy

- The most optical configurations for complete visualization of all joints
- Direct View and Videoarthroscopes delivering the versatility you need
- Excellent illumination—optimal viewing angle—sharp images—deep focus



Direct view Arthroscopes

Rapid exchange of Arthroscopes
Common interface for all camera systems
Cost effective long term ownership



VideoArthroscopes

Single handed operation
Easy reach to the light post
Fewer lenses, higher clarity



Couplers

Accepts all direct view endoscopes
Allows multi-specialty standardization
Autoclavable and non-autoclavable versions

Ordering Information

Autoclavable VideoArthroscopes for Smith & Nephew Autoclavable 470H Camera Heads

7210716	4.0 mm x 30° Autoclavable VideoArthroscope	160 mm working length
7209401	4.0 mm x 70° Autoclavable VideoArthroscope	160 mm working length
7209402	2.7 mm x 30° Autoclavable VideoArthroscope	120 mm working length
7205929	2.7 mm x 30° Autoclavable VideoArthroscope	67 mm working length

Direct View Autoclavable Arthroscopes

3894	4.0 mm x 30° Autoclavable Arthroscope	160 mm working length
3895	4.0 mm x 70° Autoclavable Arthroscope	160 mm working length
4662	4.0 mm x 0° Autoclavable Arthroscope	160 mm working length
7205681	2.7 mm x 30° Autoclavable Arthroscope	120 mm working length
7205682	2.7 mm x 30° Autoclavable Arthroscope	67 mm working length

VideoArthroscopes—Non autoclavable—for all C-mount Camera Heads

7205922	4.0 mm x 30° VideoArthroscope	160 mm working length
3628	4.0 mm x 70° VideoArthroscope	160 mm working length
3632	2.7 mm x 30° VideoArthroscope	120 mm working length
3768	2.7 mm x 30° VideoArthroscope	67 mm working length
3797	2.7 mm x 70° VideoArthroscope	70 mm working length
3863	1.9 mm x 30° VideoArthroscope	65 mm working length

Direct View Arthroscopes—Non autoclavable

3775	4.0 mm x 0° Arthroscope	160 mm working length
3623	4.0 mm x 70° Arthroscope	160 mm working length
3765	4.0 mm x 30° Arthroscope	160 mm working length
3625	2.7 mm x 30° Arthroscope	120 mm working length
3947	2.7 mm x 30° Arthroscope	120 mm working length
3948	2.7 mm x 70° Arthroscope	70 mm working length
4184	1.9 mm x 30° Arthroscope	65 mm working length

Couplers

7210192	Autoclavable coupler, for use with Smith & Nephew Autoclavable 470H camera heads and Direct View Arthroscopes
7204613	Camera coupler, C-mount, 25 mm length (black focus ring)
7204823	Camera coupler, C-mount, 30 mm length (gold focus ring)
7204614	Camera coupler, C-mount, 35 mm length (blue focus ring)
7205461	Zoom coupler, C-mount, 25-45 mm length

Endoscopy

Smith & Nephew, Inc.
Andover, MA 01810
USA

T +1 978 749 1000
F +1 978 749 1108
www.smith-nephew.com

©2003, 2006 Smith & Nephew. All rights reserved.
Printed in USA 12/06 1061191 Rev.B