



 **smith&nephew**
ELITE[◇] PASS
Premium Arthroscopic
Suture Shuttle

Trusted partner, simplified procedure



The ELITE[◇] PASS Premium Arthroscopic Suture Shuttle offers arthroscopic shoulder surgeons a flexible, single-handed method of passing suture through tissue without the need for additional instruments. The reusable handle features the trusted ACUFEX[◇] instrument design, while the disposable needle is composed of highly flexible alloy. The ELITE PASS design also gives surgeons the ability to pass suture through a single access portal, thus simplifying the procedure and decreasing costs for providers and patients.



Ordering information

- 7211020 ELITE PASS Premium Arthroscopic Suture Shuttle with Ratchet
- 7210693 ELITE PASS Shuttle Needle – box of 5