



A lasting solution in the repair of cartilaginous defects

The MOSAICPLASTY procedure is easy to learn and use with its systematic, color-coded approach

Instruments permit the flexibility of an arthroscopic or open technique in the knee and ankle

Most components are reusable or multiple-use, and the procedure can be performed in a single operation



A breakthrough technique, MOSAICPLASTY® Autogenous Osteochondral Grafting (AOCG) is the result of more than a decade of animal and clinical experience, and refinements in surgical technique and instrumentation. In clinical use since 1992, the MOSAICPLASTY procedure has been shown to be an efficacious, reproducible, and cost-effective means of restoring chondral defects. This procedure is indicated for the treatment of focal, full thickness, cartilage lesions of up to 9 cm in the knee and ankle.

Ordering Information

REF Description

Systems

7205532	MOSAICPLASTY Complete System
7205605	MOSAICPLASTY Precision System

Graft Harvesting

Tubular Chisels (five to ten patient uses)	
7207099†	2.7 mm Tubular Chisel
7207098†	3.5 mm Tubular Chisel
7207097†	4.5 mm Tubular Chisel
7205493	6.5 mm Tubular Chisel
7205494	8.5 mm Tubular Chisel

Chisel Guards

7207208†	2.7 mm Chisel Guard
7207209†	3.5 mm Chisel Guard
7207210†	4.5 mm Chisel Guard
7205499	6.5 mm Chisel Guard
7205500	8.5 mm Chisel Guard

Trephines (optional—for hard bone)

7207089	2.7 mm Trephine
7207088	3.5 mm Trephine
7207087	4.5 mm Trephine
7205497	6.5 mm Trephine
7205498	8.5 mm Trephine
7207103	Trephine Adaptor 2.7, 3.5, 4.5 mm
7205515	Trephine Adaptor 6.5, 8.5 mm

Harvesting Tamps

7207107†	2.7 mm Harvesting Tamp
7207106†	3.5 mm Harvesting Tamp
7207105†	4.5 mm Harvesting Tamp
7205495	6.5 mm Harvesting Tamp
7205496	8.5 mm Harvesting Tamp

Graft Placement

Drill Guides (five to ten patient uses)	
7207214†	2.7 mm Drill Guide
7207215†	3.5 mm Drill Guide
7207216†	4.5 mm Drill Guide
7205501	6.5 mm Drill Guide
7205502	8.5 mm Drill Guide

Drill Bits (single use)

7205508	2.7 mm Drill Bit
7207212	3.5 mm Drill Bit
7207213	4.5 mm Drill Bit
7205503	6.5 mm Drill Bit
7205504	8.5 mm Drill Bit

Dilators

7205509†	2.7 mm Dilator
7205510†	3.5 mm Dilator
7205511†	4.5 mm Dilator
7205512	6.5 mm Dilator
7205513	8.5 mm Dilator

Delivery Tamps

7207205†	2.7 mm Delivery Tamp
7207206†	3.5 mm Delivery Tamp
7207207†	4.5 mm Delivery Tamp
7205505	6.5 mm Delivery Tamp
7205506	8.5 mm Delivery Tamp

Accessories

7205507	MOSAICPLASTY Complete—Sterilization Tray Only
7205604	MOSAICPLASTY Precision—Sterilization Tray Only

MOSAICPLASTY DP Disposable System

Requires the use of the graft placement drill guides, dilators, and delivery tamps. Each set includes a disposable chisel, tamp, and drill bit.

7209234	3.5 Disposable Harvesting Set
7209235	4.5 Disposable Harvesting Set
7209236	6.5 Disposable Harvesting Set
7209237	8.5 Disposable Harvesting Set

† MOSAICPLASTY Precision System Components



Defect site after filling with ten 4.5 mm plugs

Five-year follow-up arthroscopy of previous open MOSAICPLASTY procedure



Color-coded by size

Endoscopy
Smith & Nephew, Inc.
Andover, MA 01810
USA

www.smith-nephew.com
978 749 1000
978 749 1108 Fax
800 343 5717 US Customer Service

®Trademark of Smith & Nephew.
©2005 Smith & Nephew. All rights reserved.
Printed in USA 03/05 10600007 Rev. A