

A difference you can see

Great visibility, superior absorbency<sup>1,2</sup>

 smith&nephew

**OPSITE<sup>®</sup>**  
**POST-OP VISIBLE**

Waterproof, Bacteria-proof  
Dressing with See-through  
Absorbent Pad



No hidden problems,  
just great results

Great visibility, superior absorbency<sup>1,2</sup>

 **smith&nephew**

**OPSITE<sup>®</sup>**  
**POST-OP VISIBLE**

Waterproof, Bacteria-proof  
Dressing with See-through  
Absorbent Pad




A difference you can see

Great visibility, superior absorbency<sup>1,2</sup>

 **smith&nephew**

**OPSITE<sup>®</sup>**  
**POST-OP VISIBLE**

Waterproof, Bacteria-proof  
Dressing with See-through  
Absorbent Pad



# An absorbent window for post-op sites

Great visibility, superior absorbency<sup>1,2</sup>

 **smith&nephew**

**OPSITE<sup>◇</sup>**  
**POST-OP VISIBLE**

Waterproof, Bacteria-proof  
Dressing with See-through  
Absorbent Pad



# OPSITE<sup>◇</sup> Post-Op Visible, a difference you can see

## No hidden problems

The new OPSITE Post-Op Visible dressing, with its unique lattice structure foam pad and transparent top film, has been specifically designed to give you the ability carry out regular assessments of the wound without the need to lift or remove the dressing. This allows you to react early to signs that could potentially delay the healing process.

## Just great results

Every feature of OPSITE Post-Op Visible provides tangible benefits and is designed with your patients' well being, comfort and recovery in mind. The dressing creates an optimum environment for healing by controlling moisture balance, absorbing wound discharge and preventing the build up of fluid under the dressing. This minimises both the potential for maceration and the number of dressing changes needed. Encouraging healing and reducing the chance of infection can lead to shorter recovery times and lower treatment costs.

### Bacterial barrier

Provides a barrier against bacteria, including MRSA, to reduce the risk of infection

### Lattice foam structure

Allows you to monitor progress of the wound as often as you like without unnecessary dressing changes for your patients. This reduces disturbance to the wound bed, minimises the risk of infection and helps to provide a cost effective solution

### Comfortable and conformable

The unique dressing design and low-allergy adhesive ensure it's comfortable to wear and easy to apply and remove

### Waterproof film

Patient can shower with dressing in situ

### High Moisture Vapour Transmission Rate (MVTR)<sup>2</sup>

The film has a unique molecular structure which gives a high MVTR. This allows unwanted moisture to transpire and helps prevent infection

### Highly absorbent lattice pad

Allows lateral spread across the lattice pad to maximize absorbency without impairing visibility, reducing the risk of skin maceration and minimising the number of dressing changes for cost effective care

### Ordering information

OPSITE Post-Op Visible is designed for low to moderate exuding post-operative wounds, lacerations, cuts, abrasions and minor burns.

Each OPSITE Post-Op Visible dressing is individually wrapped and sterile.

With OPSITE Post-Op Visible there are no hidden problems - just great results

**Wound Management**  
Smith & Nephew  
Medical Ltd  
101 Hessele Road  
HULL HU3 2BN

[www.smith-nephew.com/wound](http://www.smith-nephew.com/wound)

#### References

- 1 Hammond, V; DS/07/209/P1
- 2 Tompkins, L; DS/07/115/R1

T 44 (0) 1482 225181  
F 44 (0) 1482 673106

<sup>2</sup>Trademarks of Smith & Nephew  
© Smith & Nephew Oct 2007  
10680/OPOV/SINGLE SHEET/OCT07

### OPSITE Post-Op Visible Ordering Codes

Code	Size	Carton
66800137	15cm x 10cm	20
66800138	20cm x 10cm	20
66800139	25cm x 10cm	20