



**PERI-LOC<sup>◇</sup>**

Periarticular Locked  
Plating System

## Target less invasive procedures


The PERI-LOC Locked Plating Targeter enables precise screw placement for a more minimally invasive procedure. The Targeter provides accurate alignment with the plate, so that both locking and non-locking screws can be inserted percutaneously.



# PERI-LOC<sup>®</sup> Targeter System for 3.5mm Medial Distal Tibia Locking Plate

## The PERI-LOC System offers you these advantages

- Surgeon-friendly instrumentation features standardized drill bits, drill guides and screwdrivers to save time
- Precise, anatomical contouring
- Option of locking or non-locking screws in every hole



## Catalog Information

### PERI-LOC 3.5mm Targeter Pin Wire Drill Set #7181-0107

Cat. No.	Description	Qty
7117-3300	Targeter K-Wire 1.6mm X 260mm	6
7117-3406	Targeter 2.7mm PF Pin 40mm	2
7117-3418	Targeter 2.7mm Drill Bit	2
7117-3438	Targeter 2.7mm PF Pin 18mm	2

### PERI-LOC 3.5mm Medial Distal Tibia Targeter Set #7181-0103

Cat. No.	Description	Qty
7112-9401	Small Outer Case - 2.4"	1
7112-9402	Lid For Outer Cases	1
7117-0325	3.5mm D Tib Tgr Tray	1
7117-3410	Targeter 4.7mm Hexdriver	2
7117-3419	3.5mm Targeter Screw Guide	4
7117-3420	3.5mm Targeter 2.7mm Drill Guide	4
7117-3421	Targeter 1.6mm K-Wire Guide	2
7117-3422	3.5mm Targeter Trocar	2
7117-3486	3.5mm Targeter Self-Retaining Hexdriver	2
7117-3429	3.5mm Targeter Medial Distal Handle, Left	1
7117-3430	3.5mm Targeter Medial Distal Handle, Right	1
7117-3431	3.5mm Targeter Medial Distal Base, Left	1
7117-3432	3.5mm Targeter Medial Distal Base, Right	1
7117-3437	Targeter 3.5mm Base Plug	10
7117-3445	3.5mm Targeter D Tib Locking Screw	4
7117-3547	Large Screwdriver Handle	1

Orthopaedics - Smith & Nephew, Inc., 1450 Brooks Road, Memphis, TN 38116, USA  
 Telephone: 1-901-396-2121, Information: 1-800-821-5700, Orders/inquiries: 1-800-238-7538  
 www.smith-nephew.com