

# TRIGEN SURESHOT



TM

**OCHIRE DISTALA  
SIGURA**

**RAPIDA  
FARA IRADIERE**



# SISTEME MECANICE DE OCHIRE DISTALA




# OCHIREA DISTALA LA TIJ ELE ZAVORATE

- **Dureaza mult**
- **Este destul de complicata**
- **Rezultate nu intotdeauna bune**
- **Necesita experienta**
- **Expunere la radiatii**
- **Diferita in functie de tehnica**



# SOLUTIA : sistemul **TRIGEN™** **SURESHOT™**

100% sigur  
0%  
iradiere



# TRIGEN™ SURESHOT™

## Siguranta


- Imagine virtuala
- Procedura rapida

## Comoditate

- Reduce dependenta de
- Tehnica simpla, intuitiva


## Precizie

- Pozitia pacientului nu conteaza
- Feedback vizual in timp real



# CE ESTE SI CUM FUNCTIONEAZA ?

- **Tehnologie de localizare prin camp electromagnetic**
- **Senzor pe conductorul introdus in tija**
- **Imaginea virtuala a capatului distal al tije este proiectata pe ecran**
- **Asigura feedback vizual in timp real fara ajutorul razelor X**
- **Tehnica “cercurilor perfecte” - cercuri concentrice (verde si rosu)**



# INTERFATA

- **Nesterila**
- **Touch screen**
- **Nu necesita tastatura**
- **Trecere automata din plan frontal in plan sagital**
- **Platforma deschisa**
- **Timp de pornire = 1.5 minute**



# INSTRUMENTARUL SURESHOT™





# CONTROL COMPLET LA ZAVORAREA DISTALA


TRIGEN SURESHOT system: At your command

- **primul si singurul sistem cu imagine virtuala in timp real**
- **nu este o tehnica noua, ci o noua tehnologie**
- **control complet reduce iradierea scade timpul operator**

Computer-based calibrated software provides perfect circle targeting

Continuous visual real-time feedback of drill position ensures correct direction and angle

All under your direct control; requires no fluoroscopy during distal locking



# CE TIJE SE POT FOLOSI CU ACEST SISTEM ?

- **TRIGEN™ META-NAIL™ Tibia**
- **TRIGEN META-NAIL Retrograde**
- **TRIGEN TAN™**
- **TRIGEN Adolescent TAN**
- **TRIGEN FAN**
- **TRIGEN™ Humeral Nail**

