



 **smith&nephew**
TWINFIX[®] PK FT
Suture Anchor

Fully threaded anchor can be used for medial-row, lateral-row, multi-row and bridging techniques



Performance, Consistency and Simplicity.

The legacy continues with the new TWINFIX[®] PK FT Suture Anchor. Building on the clinically proven TWINFIX design, this fully threaded suture anchor made of PEEK-OPTIMA[®] polymer from Invibio[®] offers secure cortical fixation as well as sub-cortical bridging, delivering confidence in the weakest of bone conditions. The anchor boasts a new double helix thread pattern, deep eyelet placement, and smooth sliding of the ULTRABRAID[®] sutures.

Benefits

Exceptional pullout strength*

Cortical fixation and internal suture capture provide strength in the weakest bone conditions


Proven ULTRABRAID Suture

Ultra high strength suture reduces the chance of suture breakage; suture differentiation with three colors.

Robust inserter handle

TWINFIX[®] handle ensures smooth suture deployment; external posts for suture protection and control; and a secure anchor – inserter interface

**Data on file at Smith & Nephew*



Ordering Information

Two ULTRABRAID Suture Configurations (COBRAID[®] blue and COBRAID black suture)

72201576	TWINFIX [®] PK FT 5.5 mm Anchor, 2 Sutures
72201578	TWINFIX PK FT 5.5 mm Anchor, 2 Sutures with Disposable Awl
72201681	TWINFIX PK FT 6.5 mm Anchor, 2 Sutures
72201683	TWINFIX PK FT 6.5 mm Anchor, 2 Sutures with Disposable Awl

Three ULTRABRAID Suture Configurations (COBRAID blue, COBRAID black and white suture)

72201577	TWINFIX PK FT 5.5 mm Anchor, 3 Sutures
72201579	TWINFIX PK FT 5.5 mm Anchor, 3 Sutures with Disposable Awl
72201682	TWINFIX PK FT 6.5 mm Anchor, 3 Sutures
72201684	TWINFIX PK FT 6.5 mm Anchor, 3 Sutures with Disposable Awl