

 smith&nephew

**TWINFIX<sup>®</sup>**

Suture Anchors with  
ULTRABRAID<sup>®</sup> Suture

Unparalleled strength,  
superior handling



The confidence of a strong, secure shoulder repair begins with ULTRABRAID<sup>®</sup> Suture. This non-absorbable suture is composed of a special ultra high molecular weight (UHMW) polyethylene fiber and features a unique braid configuration. As a result, ULTRABRAID Suture offers special advantages over traditional polyester suture, including higher knot breaking strength, increased lubricity, and a stronger resistance to fraying<sup>1</sup>. ULTRABRAID Suture is preloaded on all TWINFIX<sup>®</sup> Suture Anchors that have #2 sutures.

In addition to the ULTRABRAID Suture, TWINFIX Suture Anchors have a patented Hi-Lo thread geometry to help attain consistent failure loads in all types of bone quality.

TWINFIX AB Suture Anchors also provide a tapered body that is mated to the dilator to reduce insertion torque, thereby helping to minimize anchor breakage. Suture eyelets are located deeper in the anchor body to reduce suture eyelet failure under cyclic loading.

TWINFIX Ti Suture Anchors provide secure fixation with a wide range of sizes to match various applications and bone qualities. Anchor cores are matched to the drill size to facilitate insertion while maximizing fixation strength.

For arthroscopic or mini-open shoulder procedures, trust the strength and security of TWINFIX Suture Anchors with ULTRABRAID Suture.



Titanium and Absorbable Anchors

## Ordering Information

### Titanium Anchor

7210595	TWINFIX Ti 2.8 mm Suture Anchor with one #2 ULTRABRAID Suture
7210597	TWINFIX Ti 3.5 mm Suture Anchor with two #2 ULTRABRAID Sutures
7210598	TWINFIX Ti 5.0 mm Suture Anchor with two #2 ULTRABRAID Sutures
7210599	TWINFIX Ti 6.5 mm Suture Anchor with two #2 ULTRABRAID Sutures
7210707	TWINFIX Ti 3.5 mm Suture Anchor with two #2 ULTRABRAID Sutures and needles
7210708	TWINFIX Ti 5.0 mm Suture Anchor with two #2 ULTRABRAID Sutures and needles
7210709	TWINFIX Ti 6.5 mm Suture Anchor with two #2 ULTRABRAID Sutures and needles



Titanium Anchor

### Absorbable Anchor

7210701	TWINFIX AB 5.0 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures
7210702	TWINFIX AB 5.0 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures, packaged with threaded dilator and guide
7210703	TWINFIX AB 6.5 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures
7210704	TWINFIX AB 6.5 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures, packaged with threaded dilator
7210710	TWINFIX AB 5.0 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures and needles
7210711	TWINFIX AB 5.0 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures, needles and dilator
7210712	TWINFIX AB 6.5 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures and needles
7210713	TWINFIX AB 6.5 mm Suture Anchor, PLLA, with two #2 ULTRABRAID Sutures, needles and dilator

### Accessories

7209500	Drill Guide, Spiked Tip – 6" x 3.5 mm
7209115	Drill Guide, Spiked Tip – 6" x 5.0 mm
7210269	Drill Guide, Spiked Tip – 6" x 6.5 mm
7209506	1.8 mm Drill for 2.8 mm Anchor
7209502	2.5 mm Drill for 3.5 and 5.0 mm Anchor
7210215	3.5 mm Drill for 6.5 TWINFIX Ti Anchor
7209298	5.0 mm Threaded Dilator
7210383	6.5 mm Threaded Dilator

<sup>1</sup>Data on file at Smith & Nephew