


 **smith&nephew**
LEUKOSTRIP[®]
Adhesive Hypoallergenic
Wound Closure Strips

Flexible, adhesive, hypoallergenic wound
closure strips



LEUKOSTRIP[◇], is flexible, adhesive, hypoallergenic wound closure strips

LEUKOSTRIP wound closure strips are made of 100% polyamide material coated with a hypoallergenic adhesive which allows almost pain free closure of the wound. By replacing skin clips with LEUKOSTRIP on the third post-operative day, scar formation at the puncture sites is minimised. LEUKOSTRIP can also be combined with subcutaneous or intracutaneous stitches as well as being used to supplement clip or individual over-and-over sutures.

Advantages

Flexible

The polyamide material is elastic. After closing the wound it returns to its original length.

Helps protect the wound against tensile and shear forces

LEUKOSTRIP does not traumatise the edges of the wound and allows an almost pain free closure of the wound as well as minimising scar formation.

Permeable

LEUKOSTRIP is permeable to water vapour and air and thus reduces the risk of maceration.

Hypoallergenic

The adhesive is hypoallergenic, reducing the risk of allergic reactions.

Reliable adhesive strength

The dressing adheres, reliably and securely.


Indications

Primary and secondary closure of wounds

Fixing skin grafts

Replacing skin sutures following subcutaneous sutures

Support and relief for intracutaneous sutures, clip stitching or individual over-and-over sutures.



Ordering Information

Suture material boxes

Code	Size	Contents Strips/Bag	Contents per Bag Unit	Contents per Packaging Unit
66202876	38mm x 4.0mm	8	50	4 x 50
66002878	76mm x 6.4mm	3	50	4 x 50
66002879	102mm x 6.4mm	5	50	4 x 50
66002880	102mm x 13.0mm	6	50	4 x 50
66002881	102mm x 26.0mm	4	25	4 x 25

Cartons

Code	Size	Contents Strips/Bag	Contents per Bag Unit	Contents per Packaging Unit
66072954	38mm x 4.0mm	8	10	20 x 10
66002952	76mm x 6.4mm	3	10	20 x 10
66072955	102mm x 6.4mm	5	10	20 x 10
66002953	102mm x 13.0mm	6	10	20 x 10

Wound Management www.smith-nephew.com
Smith & Nephew
Medical Ltd
101 Hessle Road
Hull HU3 2BN
England

T 44 (0) 1482 225181
F 44 (0) 1482 328326

[◇]Trademarks of Smith & Nephew
6853/LEUKOSTRIP/SS/GLOBAL/0404